

# Good Vibrations™

## Ultrasonic Cleaner Solution



*Clean and protect  
with Good Vibrations.*

Good Vibrations is a super concentrated instrument cleaning solution formulated for use in ultrasonic cleaning devices. It provides rapid and thorough removal of a wide range of soils while protecting delicate instruments.

### Enhances Ultrasonic Cleaning Action

- Non-foaming formula best for ultrasonic cleaners.
- Penetrating agents improve cleaning action in small crevices and hinges.
- Odor suppressant neutralizes organic and chemical odors.

### Protects Delicate Instruments

- Team of anti-corrosive agents protects metal surfaces.
- Best results for carbon steel items including burs and curets.
- Safe cleaning for jointed orthodontic pliers.
- Enhances instrument appearance and shine.

### Biodegradable and Contains No Phosphates.

### Economical Super Concentrate – One Ounce Makes One Gallon of Cleaning Solution

### Meet CDC Guidelines for Instrument Decontamination

- Effective cleaning is essential prior to sterilization.
- Effective ultrasonic cleaning reduces risk of injury during manual scrubbing of instruments.



Good Vibrations Description	Case Qty	Order #
1 oz. Unit Dose - (24/box)	6	GVUS24
1 Gallon with Pump	4	GVUS128
55 Gallon Drum	1	GVUS550

## Good Vibrations™

### Ultrasonic Cleaner Solution

#### Description:

Good *Vibrations* is a super concentrated instrument cleaning solution formulated for use in ultrasonic cleaning devices. The product is available in one ounce unit dose packets and in gallons with metered pump. It may be diluted with ordinary tap water.

#### Intended Use:

Good *Vibrations* is intended for cleaning instruments by aiding the cavitation process of ultrasonic cleaners. It is designed to minimize manual handling/scrubbing of contaminated items and is especially useful for cleaning instruments with serrated edges, joints and hinges.

#### Product Specifications:

Dilution Ratios: 1 oz. or 1 unit dose per gallon.  
Removal Time: 5 to 12 minutes depending on volume and type of soil.  
pH Concentrate: 11.0  
Use Temperature: Room temperature, 68°F (20°C).  
Odor: Mint  
Shelf Life: 3 years from date of manufacture.  
Disposal: Dispose according to all local, state and federal regulations. This product contains no phosphates.

#### Product Capabilities:

Good *Vibrations* rapidly removes organic and inorganic soils. The detergent formula lifts and disperses soils during the ultrasonic cleaning process.

A powerful anti-corrosive system protects delicate instruments, including carbon steel and anodized aluminum. The anti-corrosive agents are especially effective to protect box locks in jointed items and hinges.

#### Use Precautions:

Good *Vibrations* is an eye irritant and mild skin irritant. Follow CDC Guidelines for protective attire during instrument decontamination: Safety goggles and/or face shield, mask (for biological hazards) and water-resistant gloves. Solution will fade when exposed to sunlight. Store in a cool, dry area away from direct sunlight.

Due to the corrosive effects of water, it is not recommended to soak metal items for prolonged periods (over two hours or overnight).

#### Accessory Products:

Other products available from Certol include:



**ProSpray™**  
Ready-to-Use Surface  
Disinfectant/Cleaner



**ProSpray™ wipes**  
Surface Disinfectant/  
Cleaner Towelettes



**ProEZ foam™**  
Ready-to-Use Foaming  
Enzymatic Spray



**ProLube RTU™**  
Instrument Lubricant  
Ready-to-Use