

SUPER CONCENTRATED DETERGENT LINE

Certol[®]
International

ProEZ SCD™

Super Concentrated Quad Enzyme Automatic Washer Detergent

Non-Foaming • Corrosion Inhibitors • Extended High Temperature Performance



Super concentrated enzyme cleaning power for all surgical soils

ProEZ SCD is part of the:

SUPER CONCENTRATED DETERGENT LINE

SUPER SYSTEM



Enzymatic



Neutral



Alkaline



Lubricant

Super Concentrated:

- Clinically tested for effective cleaning at ultra low dose settings.
- Reduces weight of containers and storage space.
- Delivers critical cleaning performance WITH economy.

Stabilized Four Enzyme Formula:

- Specific enzymes cleans proteins/blood, fats, starches and sticky components of biofilm.
- Improves cleaning results for tough orthopedic cases.
- Effective at expanded temperature range for both enzymatic AND regular wash cycles up to 150°F.
- Preservatives maintain enzyme effectiveness over typical storage periods.

Enhanced Cleaning Action in Automated Equipment:

- Non-foaming for easy rinsing.
- Promotes ideal pump pressure and spin action for best cleaning action.

Broad Material Compatibility:

- Extensive testing demonstrates compatibility with stainless steel, titanium, gold, brass, copper, aluminum and anodized aluminum.
- Neutral pH and lab tested anti-corrosive system.



ProEZ SCD Description	Case Qty	Order#
2.5 Gallon Drum	1	EZSCD2.5

ProEZ SCD™

Super Concentrated Quad Enzyme Automatic Washer Detergent

Description:

ProEZ SCD Super Concentrated Quad Enzyme Automatic Washer Detergent is non-foaming and neutral pH. The multi-enzyme formula quickly dissolves all common surgical soils including proteins, blood, fats and starches. It is super concentrated to reduce package weight and storage space.

Product Capabilities:

ProEZ SCD contains protease, amylase, lipase and cellulase enzymes to rapidly break up a wide range of surgical soils including blood, fats, orthopedic mucous and starches. It is non-corrosive and safe to use on instruments and devices indicated for cleaning in automated washing equipment.

Product Specifications:

Dilution Ratios: 1/40 to 1/5 oz. per gallon
(0.2 ml to 1.6 ml per liter).
pH: Neutral.
Temperature of Use: 90°F to 150°F (32°C to 66°C).
Odor: Non-specific, fresh.
Color: Natural clear yellow.
Shelf Life: 2 years from date of manufacture.

Intended Use:

ProEZ SCD Super Concentrated Quad Enzyme Automatic Washer Detergent is intended for use when cleaning instruments and appropriate devices. This product may be used in ultrasonic cleaning devices, AERs (automated endoscope reprocessors) and automated washers at low and high impingement in prewash, enzyme and main wash cycles. It is compatible with a wide range of materials including regular and anodized aluminum, gold, copper, brass, titanium, plastics, glass optics and stainless steel. All applications should be according to manufacturer's instructions.

Accessory Products:

Other products available from Certol include:



Certol Family of Super Concentrated Detergents

ProEZ SCD™

- QUAD Enzymatic

ProWash SCD™

- Neutral

ProClense plus SCD™

- Alkaline

ProLube SCD™

- Lubricant



ProSpray™ wipes

- Surface Disinfectant/
• Cleaner Towelettes



ProEZ foam™

- Ready-to-Use Foaming
• Enzymatic Spray