

STEP 1

Preclean at Point-of-Use

- Dried soils on instruments?
- Spills during transport with presoak pans?
- Labor and Delivery creating overnight cleaning backlog?



Solution:

ProEZ foam™
 Spray ProEZ foam and walk away!
 Let enzymes do the work for you.

1

STEP 2

Enzyme Soak/Scope Cleaning

- Need a powerful, yet safe solution?

Solution Choices:

- ProEZ AW Quad™**: Premium four enzyme
- or
- ProEZ 2™**: Dual enzyme
- or
- ProEZ 1™**: Economical single enzyme.



2

The Certol® System
 Customize Your System
 According to Your
 Challenges

STEP 3

Clean and Rinse

- Need an instrument detergent safe for all metals?

Solution Choices:

- ProWash™**: Premium detergent and cart washer
- or
- ProClense plus™**: Economical choice for a detergent with advanced metal conditioning agents.



3

STEP 4

Lubricate and Protect

- Joints stiff? Instrument joints too?
- Infection control wants that "milk bath" out of action?
- Looked at your current lubricant label lately? How do you shake 15 gallons?

Solution:

ProLube™ Concentrate or ProLube™ RTU
The Clear Choice™
 Clearly better, oil-free and non-separating.



4

STEP 5

Sterilize

REMEMBER!
 Sterile is not sterile if clean is not clean. After thorough cleaning with Certol products, proceed to your regular inspection and sterilization of your surgical instruments.



5