



## Technical Bulletin: Use of ProSpray disinfectants on blood pressure cuffs

- ProSpray disinfectants are an intermediate level (inactivates tuberculosis var. bovis) surface disinfectant. They are water-based and meet manufacturer's recommendations for use of an alcohol free disinfectant.
- ProSpray disinfectants may be used on a wide variety of equipment including blood pressure cuffs.

ProSpray water-based intermediate level disinfectants are available in the following formats and corresponding order codes:

- ProSpray ready to use liquids: PSC240 (24 oz. spray bottle); PSC128 (gallons); PSC050 (5 gallon drum).
- ProSpray wipes disinfectant towelettes: PSW (soft pack); PSWC (standard canister); PSWCT (tall canister).
- ProSpray C60 concentrate: PSC60/48 (unit dose packets); PSC60/32 (quart tip & pour).